

Cryo Chambers & Beds

Dominium Group d.o.o.
www.dominiumgroup.rs
office@dominiumgroup.rs
+381 60 46 222 22



DOMINIUM
WELLNESS

Introduction



Cryotherapy — advanced cold therapy that reduces inflammation, relieves pain, and primes your body for what's next. Cryotherapy chambers capable of reaching temperatures as low as -190°C are at the forefront of cutting-edge therapies for physical recovery, pain management, and performance enhancement. This advanced treatment harnesses extreme cold to stimulate the body's natural healing processes, offering a wealth of benefits for both patients managing health conditions and athletes striving to perform at their peak.

Benefits of Cryotherapy for Athletes



1.Speeds Muscle Recovery

Cryotherapy is a go-to treatment for elite athletes due to its ability to significantly reduce muscle soreness and stiffness after intense training or competition. The therapy flushes out lactic acid buildup, enhancing muscle recovery and performance.



2.Enhances Athletic Performance

Regular cryotherapy sessions can improve endurance and strength by optimising recovery time and maintaining peak physical condition. Athletes experience reduced downtime, allowing them to train more efficiently.



3.Alleviates Injuries and Prevents Overuse

Common sports injuries like sprains, strains, and inflammation benefit from cryotherapy's anti-inflammatory effects, enabling athletes to recover faster and return to their sport.



4.Boosts Energy and Focus

The cold environment stimulates the release of adrenaline and other natural chemicals, enhancing energy, focus, and mental clarity—critical factors for peak



DOMINIUM
WELLNESS

DWKS80E

Cryo Chamber Specs

NO.	item	Liquid Nitrogen Chamber
1	Temperature Range	-80°C~Room Temp
2	Cooling Time (to -80°C)	20min
3	Ext. Dimensions (W* D* H)	2000*900*2350mm
4	Int. Dimensions (W* D* H)	900*700*2050mm
5	Maximum User Weight	0-250KG
6	Capacity	1-2 People
7	Power specification	220V/50Hz
8	Installed Power	8KW
9	Cooling Source	Cascade Compressor Unit
10	Cooling Speed	Reaches -80°C within 20 minutes
11	Install	All-in-One, Plug-and-Play, with no installation required



DOMINIUM
WELLNESS

DWKK140A

Cryo Chamber Specs

Ultra-Low Cryo Chamber Specs		
NO.	Item	DWKK140A
1	Temperature Range	-140 °C~ Room Temperature
2	Dimensions (W ×D×H)	1000 * 1400 * 2200mm
3	User Height Range	150-220cm
4	Maximum User Weight	0-150kg
5	Device Weight	about 300kg
6	Nitrogen Use per Session	3-5KG
7	Power Supply Requirement	120V/60HZ Or Customizable
8	Installed Power	≤ 2KW
9	Cooling Source	Liquid Nitrogen

Additional Features

Ultra-Low Cryo Chamber – Main Components		
NO.	Item	DWKK140A
1	Dual-Screen Display	√
2	Electric Lifting Platform	√
3	Thermal Therapy System	/
4	Control System	√
5	Forced Convection Cold Air System	√
6	Indoor Nitrogen Exhaust System	/
7	Ambient Lighting	√
8	Intelligent Cold Air Control System	/
9	Warranty	1 year